










Produktnamn: **biovalvesafe®**

Art.nr	Yttre Ø	Inre Ø	G	Längd	Flödes hastighet (passivt)	Flödes hastighet Tryck 350 PSI	Färgkod
106.082	0.9 mm	0.6 mm	22	25 mm	36 ml/min	7 ml/s	
106.102	1.1 mm	0.8 mm	20	32 mm	60 ml/min	11,9 ml/s	
106.1232	1.3 mm	1.0 mm	18	33 mm	97ml/min	18,6 ml/s	
106.122	1.3 mm	1.0 mm	18	45 mm	90 ml/min	17,1 ml/s	
106.142	1.5 mm	1.2 mm	17	45 mm	125 ml/min	23 ml/s	
106.172	1.7 mm	1.4 mm	16	45 mm	180 ml/min	>23,3 ml/s	
106.212	2.1 mm	1.6 mm	14	45 mm	270 ml/min	>23,3 ml/s	

Avdelningsförpackning: 50 st

Transportförpackning: 500 st

Material: Polypropylene (PP), Polyetylene (PE), PTFE, och rostfritt stål. PVC-fri och latexfri

Övrig information: Stickskyddad PVK. CT-Godkänd 350 PSI

Produkten är CE-märkt: Ja